


# Cytiva™ Plus Cardiomyocytes

A relevant human model of  
cardiotoxicity

---

## Product Instructions

Codes:      ES-7001       $\geq 1.0 \times 10^6$  (1,000,000) viable cells  
              ES-7002       $\geq 3.5 \times 10^6$  (3,500,000) viable cells



---

**ASCENDANCE™**

esibio.com

East Coast: 200 Boston Avenue, Suite 1500 • Medford, MA 02155 • t 781-391-0205

West Coast: 1010 Atlantic Avenue, Suite 102 • Alameda, CA 94501 • t 510-522-2856

1

MAN-001, Rev.00

# Table of Contents

Page(s)

<b>1. Legal</b> .....	3
<b>2. Handling</b> .....	4
2.1. Safety warnings and precautions	
2.2. Storage	
2.3. Expiration	
2.4. Packaging	
<b>3. Introduction</b> .....	5
<b>4. Components and other materials required</b> .....	5-6
4.1. Components	
4.2. Materials to be supplied by user	
4.3. Equipment needed	
<b>5. Protocols for Multi-Electrode Arrays (MEA)</b> .....	6-13
5.1. Preparation of MEA coated plates	
5.2. Preparation of RPMI 1640/B27 medium for MEA	
5.3. Thawing of Cytiva Plus Cardiomyocytes for MEA	
5.4. Determining Post-Thaw Cell Viability of Cytiva Plus Cardiomyocytes for MEA	
5.5. Seeding Cytiva Plus Cardiomyocytes onto MEA plates	
5.6. Media change on day 4 post-thaw	
<b>6. Protocols for High Content Analysis (HCA)</b> .....	14-19
6.1. Preparation of cell culture plates for HCA	
6.2. Preparation of RPMI 1640/B27 medium for HCA	
6.3. Thawing of Cytiva Plus Cardiomyocytes for HCA	
<b>7. Protocols for Manual Patch Clamp (MPC)</b> .....	20-26
7.1 Sterilization of glass coverslips	
7.2 Preparation of 1:2 diluted Matrigel aliquots	
7.3 Preparation of 1:30 diluted Matrigel	
7.4 Preparation of Matrigel-coated coverslips	
7.5 Preparation of RPMI 1640/B27 medium for MPC	
7.6 Thawing Cytiva Plus Cardiomyocytes for MPC	
7.7 Determining Post-thaw Viability for MPC	
7.8 Seeding Cytiva Plus Cardiomyocytes onto coverslips	
7.9 Media change of day 4 post-thaw	
<b>8. Troubleshooting</b> .....	27-28
<b>9. Product Protocol Flowchart</b> .....	29

---

**ASCENDANCE**<sup>™</sup>

## Legal

Cytiva® is a trademark of General Electric Company. B27® and KnockOut® are trademarks of Gibco. Matrigel® is a trademark of Becton, Dickinson and Company. NucleoCounter®, NucleoCassette® and NC-100® are trademarks of Chemometec A/S. µClear® is a registered trademark of Greiner Bio-One GmbH.

All other third party trademarks are the property of their respective owner.

**License statements:** The following is important license information. The purchase of Cytiva Plus Cardiomyocytes includes a limited license to use the Cytiva Plus Cardiomyocytes for internal research and development, but not for any commercial purposes. Commercial purposes shall include: sale, lease, license or other transfer of the material or any material derived or produced from it; sale, lease, license or other grant of rights to use this material or any material derived or produced from it; use of this material to perform services for a fee from third parties. A license to use the Cytiva Plus Cardiomyocytes for commercial purposes is subject to a separate license agreement with GE Healthcare. IF YOU REQUIRE A COMMERCIAL LICENSE TO USE THIS MATERIAL AND DO NOT HAVE ONE, PLEASE CONTACT YOUR LOCAL ASCENDANCE REPRESENTATIVE AND RETURN THIS MATERIAL UNOPENED TO ASCENDANCE BIOTECHNOLOGY, 1010 ATLANTIC AVENUE, SUITE 102, ALAMEDA CA 94501.

Cytiva Plus Cardiomyocytes (the Products) are sold under license from GE Healthcare, Asterias Biotherapeutics Inc., and Wisconsin Alumni Research Foundation (WARF) under US patent and publication numbers: US 5,843,780, US 6,200,806, US 6,602,711, US 6,800,480, US 7,005,252, US 7,029,913, US 7,297,539, US 7,413,902, US 7,425,448, US 7,452,718, US 7,582,479, US 7,732,199, US 7,763,464, US 7,781,216, US 7,851,167, US 7,897,389 and US 2007/0010012 and equivalent patent and patent applications in other countries.

**Important note:** For *in vitro* research use only. Not to be used for any therapeutic or diagnostic applications. Not recommended or intended for diagnosis of disease in humans or animals. Do not use internally or externally in humans or animals.

All goods and services are sold subject to the terms and conditions of sale. A copy of these terms and conditions is available on request. Contact your local Ascendace representative for the most current information.

<http://www.esibio.com>

© 2017 Ascendace Biotechnology – All rights reserved.

---

**ASCENDANCE™**

## 2. Handling

Cytiva Plus Cardiomyocytes are cryopreserved in a solution containing dimethyl sulphoxide (DMSO). A Material Safety Data Sheet (MSDS) for the DMSO is available upon request by phone ((510)-522-2856) or by email at [techsupport@esibio.com](mailto:techsupport@esibio.com).

### 2.1. Safety warnings and precautions

**Warning: For research use only.** Not recommended or intended for diagnosis of disease in humans or animals. Do not use internally or externally in humans or animals. All chemicals should be considered as potentially hazardous. We therefore recommend that this product be handled only by those persons who have been trained in laboratory techniques and that it is used in accordance with the principles of good laboratory practice. Wear suitable protective clothing such as laboratory overalls, safety glasses and gloves. Care should be taken to avoid contact with skin or eyes. In the case of contact with skin or eyes wash immediately with water. See material safety data sheet(s) and/or safety statement(s) for specific advice.

### 2.2. Storage

Cytiva Plus Cardiomyocyte product codes (ES-7001) and (ES-7002) are supplied in 1 ml cryovials of  $\geq 1 \times 10^6$  (1,000,000) viable cells ( $1.6 \times 10^6$  total cells) and  $\geq 3.5 \times 10^6$  (3,500,000) viable cells ( $5 \times 10^6$  total cells) respectively. Cells are cryopreserved in 10% DMSO and 90% fetal bovine serum. Upon receipt, the frozen vial(s) of cells should immediately be removed from the outer packaging and transferred to the vapor phase of a liquid nitrogen storage unit at  $-140^\circ\text{C}$ .

### 2.3. Expiration

Cells are stable for at least 6 months when stored in liquid nitrogen (vapor phase). Cells can be maintained for at least 7 days once recovered into culture.

### 2.4. Packaging

Cytiva Plus Cardiomyocytes are provided as cryopreserved single cell suspensions in 1 mL cryovials.

### 3. Introduction

Cytiva Plus Cardiomyocytes are human cardiomyocytes derived from the NIH approved stem cell line NIH hESC-10-0061 (WA07 (H7)) and have been extensively characterized and functionally verified by flow cytometry, sub-cellular imaging and electrophysiology. Cytiva Plus Cardiomyocytes are supplied cryopreserved in a ready to use format.

### 4. Components and other materials required

#### 4.1. Components

- Cytiva Plus Cardiomyocytes
- Certificate of Analysis including product specifications
- User manual

#### 4.2. Materials to be supplied by user

<b>Product</b>	<b>Supplier/Product code</b>
• <u>15 mL Centrifuge tubes</u>	Corning 430052
• <u>50 mL Centrifuge tubes</u>	Corning 430290
• <u>Serological pipettes – 5 mL</u>	Corning 4487
• <u>Serological pipettes – 10 mL</u>	Corning 4488
• <u>Serological pipettes – 25 mL</u>	Corning 4489
• <u>500 mL Cellulose acetate filter units</u>	Corning 430769
• <u>12- or 48-well plate for Multi-Electrode Array (MEA) applications</u>	e.g. Axion BioSystems M768-GLx or M768-KAP-48
• <u>Multi-well cell plates for High Content Analysis (HCA) and Ca<sup>2+</sup> transient applications</u>	Various, e.g. 96-well plate Greiner µClear® 781091
• <u>Multi-well plates for impedance applications</u>	e.g. ACEA E-plate Cardio 96
• <u>Coverslips - electrophysiological applications*</u>	VWR 631-0149
• <u>1.5 mL tubes</u>	Corning 430290
• <u>Sterile bottle (e.g. 30 or 60 mL)</u>	NALG2019-0030/0060
• <u>RPMI 1640 + Glutamine</u>	Gibco 21875034
• <u>D-PBS</u>	Sigma D8537
• <u>B27 (50x)</u>	Gibco 17504-004
• <u>Fibronectin</u>	BD Biosciences 354008
• <u>FBS</u>	Gibco 26140-079
• <u>Matrigel®</u>	Becton Dickinson 356231
• <u>KnockOut® D-MEM (KO-DMEM)</u>	Gibco 10829-018
• <u>Sterile distilled water</u>	Fresenius Kabi 22-96-985

\* –not available in North America. Standard glass coverslips of approximately 13 mm diameter and 0.13 mm thickness are suitable.

### 4.3. Equipment needed

- Adjustable pipettes and tips
- Liquid nitrogen vapor store
- Biosafety cabinet
- Ice bucket with dry ice
- Cryovial rack
- Vacuum pump and line
- Hemocytometer or automated cell counter
- 37°C water bath
- Centrifuge

#### Application dependent:

- Patch clamp system
- Sub-cellular imaging system for High Content Analysis
- Multi-Electrode Array system
- Impedance system
- Calcium transient system

## 5. Protocols for Multi-Electrode Arrays (MEA)

### 5.1. Preparation of MEA coated plates (Time 4 hours)

#### Consumables

- FBS (at 4°C)
- Sterile distilled H<sub>2</sub>O
- Fibronectin (1 mg, human, BD Biosciences 354008) (at 4°C)
- 50 ml tube
- D-PBS

#### Equipment

- MEA plate(s) (e.g. Axion BioSystems 12-well or 48-well plates)\*
- 20 µL pipette & sterile 1-20 µL tips
- 1,000 µL pipette & sterile 1-1,000 µL tips


**\*Note:** Approximate number of MEA plates that can be used with different vial sizes of Cytiva Plus Cardiomyocytes (assuming  $6.00 \times 10^4$  cells/well; > 70% viability).

Cytiva Plus Cardiomyocyte Vial Size	Approx. Number of MEA Plates		
	6-well	12-well	48-well
ES-7001 - $1.0 \times 10^6$ viable cells/vial	2.8	1.4	0.3
ES-7002 - $3.5 \times 10^6$ viable cells/vial	9.7	4.9	1.2

## Protocol

Perform the following steps aseptically inside the biosafety cabinet (BSC). The first step in preparing an MEA plate for use is to ensure that the surface is hydrophilic. The surface of a new MEA plate is hydrophobic, and even hydrophilic MEA wells tend to become hydrophobic again during storage. A hydrophobic surface will prevent attachment and growth of the (hydrophilic) cells. FBS treatment renders the surface hydrophilic.

1. Place a 4  $\mu\text{L}$  bead of FBS over the recording electrode area of each well of the MEA plate (see Fig 1). **Note:** If a bead fails to form, ignore that recording site.
2. For a 12 well MEA plate, add 500  $\mu\text{L}$  of sterile distilled  $\text{H}_2\text{O}$  to the gaps between the wells to prevent the evaporation of the bead of FBS (see Fig 1). For a 48 well MEA plate add 200  $\mu\text{L}$  of sterile distilled  $\text{H}_2\text{O}$  to the gaps between the wells.
3. Put the lid on the MEA plate and incubate the MEA plate for 1.5 hours at room temperature.
4. Prepare a 1 mg/mL solution of fibronectin by adding 1 mL sterile distilled  $\text{H}_2\text{O}$  to the 1 mg fibronectin.
5. Take 12.5  $\mu\text{L}$  of this 1 mg/mL fibronectin solution and add it to 987.5  $\mu\text{L}$  D-PBS in a 50 mL tube for a final concentration of 12.5  $\mu\text{g}/\text{mL}$  fibronectin.
6. Aspirate the FBS bead from each well of the MEA plate (use a pipette set to dispense 8  $\mu\text{L}$ ) and **immediately** replace with a 4  $\mu\text{L}$  bead of 12.5  $\mu\text{g}/\text{mL}$  fibronectin solution over the recording electrode area.



**Figure 1.** MEA plate showing location of sterile water addition to the inter-well space to prevent evaporation of the FBS and fibronectin drop.

## 5.2. Preparation of RPMI 1640/B27 medium for MEA (Time 30 minutes)

### Consumables

- RPMI 1640 + Glutamine medium, 500 mL (at 4°C)
- B27<sup>®</sup> supplement, 10 mL (at 4°C)
- 10 mL sterile serological pipette

### Equipment

- 500 mL 0.22 µm cellulose acetate filter unit (Corning 430769)
- Vacuum line
- Pipette gun

### Protocol

During the 2-hour fibronectin incubation step, prepare medium.

**Note:** Medium may be prepared immediately before thawing Cytiva Plus Cardiomyocytes or prepared and stored at 2–8°C and used within one week of preparation.

1. Thaw frozen 10 mL B27 supplement vial(s) in a 37°C water bath for 10 minutes. Do not incubate at 37°C or expose to light for extended periods of time. Perform the following steps aseptically inside the biosafety cabinet (BSC).
2. Wipe the required number of RPMI 1640 + Glutamine medium bottle(s) and B27 Supplement vial(s) with 70% isopropanol and transfer to a BSC.
3. Place the filtration unit into the BSC.
4. Using a 10 mL sterile serological pipette, add 10 mL of B27 supplement to a 500 mL bottle of RPMI1640 + Glutamine medium. Swirl bottle several times to mix.
5. Carefully transfer the medium and supplement into the reservoir of the filter unit.
6. Place the lid on the filter unit.
7. Connect the filter unit to a vacuum source.
8. When filtration is complete, disconnect the vacuum source.
9. Detach the upper reservoir of the filtration unit.
10. Place sterile cap on the bottle portion of the filter unit.

**Note:** Store the medium at 2–8°C. Use within one week of preparation. Avoid repeated warming of the RPMI 1640/ B27 medium. Warm only the required volume of medium to complete the task.



### 5.3. Thawing Cytiva Plus Cardiomyocytes for MEA (Time 30 minutes)

#### Consumables

- 3.5×10<sup>6</sup> vial of Cytiva Plus Cardiomyocytes (ES-7001)
- RPMI 1640/B27 medium (filter sterilized)
- 50 mL tube
- 10 mL sterile serological pipette

#### Equipment

- Ice bucket with dry ice
- Cryovial rack
- 37°C water bath
- Centrifuge
- 1,000 µL pipette & sterile 1-1,000 µL tips
- Pipette gun

#### Protocol

During the 2-hour fibronectin incubation step, perform the following steps aseptically inside the biosafety cabinet (BSC).

1. Remove the cryovial from cryostorage and place onto dry ice until ready to thaw.
2. Thaw cell suspension in a 37°C water bath with gentle agitation until ice crystals just disappear.

**Note:** Take care not to immerse the whole cryovial into the water bath. Avoid extended incubation at 37°C.

3. Wipe the outside of the cryovial with 70% isopropanol and transfer to BSC.
4. Carefully transfer the cell suspension into a sterile 50 mL centrifuge tube using a 1,000 µL pipette.
5. Rinse the inside of the cryovial with 1 mL of room temperature RPMI 1640/B27 and combine with the cell suspension drop-wise with gentle mixing.
6. Slowly (over the course of 2 minutes) add 8 mL of RPMI 1640/B27 to the 50 mL centrifuge tube.
7. Centrifuge at 300 g for 5 minutes at 20°C.
8. Carefully, using a 10 mL stripette, remove 8 mLs of the supernatant. Remove a further 1 mL with a 1,000 µL pipette taking care not to disturb the cell pellet. Re-suspend the cells in the residual liquid (approx. 1 mL) with gentle agitation.

## 5.4. Determining Post-Thaw Cell Viability of Cytiva Plus Cardiomyocytes for MEA (Time 30 minutes)

During the 2-hour fibronectin incubation step, determine the viable cell number and viable cell density using preferred method of choice. We use a NucleoCounter NC-100 (ChemoMetec) cell counter, the method for which is outlined below. However, this protocol could be adapted for other commercial cell counters.

### Consumables

- RPMI 1640/B27 medium (filter sterilized)
- 4×1.5 mL tubes
- 4×NucleoCassettes™ (ChemoMetec)
- Reagent A100 (ChemoMetec)
- Reagent B (ChemoMetec)

### Equipment

- 100 µL pipette & sterile 1-200 µL tips
- 1,000 µL pipette & sterile 1-1,000 µL tips
- NucleoCounter NC-100 (ChemoMetec)

### Protocol

Perform the following steps aseptically inside the biosafety cabinet (BSC).

1. Determine viable cell number using the NucleoCounter NC-100 (ChemoMetec). To do this, transfer 40 µL of the cell suspension into a 1.5 mL tube. Add to the tube 360 µL of RPMI 1640/B27 medium (i.e. now a 1:10 dilution). Gently agitate to achieve an even cell suspension.
2. The remaining steps can be performed outside of the BSC.
3. Transfer 100 µL of 1:10 diluted cell suspension into each of three more 1.5 mL tubes.
4. Use 100 µL of 1:10 diluted cell suspension sample to calculate the number of **Non-Viable Cells/mL** (i.e. directly load the sample into a NucleoCassette).
5. Repeat step 4 for a second time, and calculate the average result.
6. Use another 1×100 µL of 1:10 diluted cell suspension sample to calculate the **Total Number of Cells/mL**. To do this, add 100 µL of Reagent A100, and then 100 µL of Reagent B, to the 100 µL of 1:10 diluted cell suspension sample. Mix by pipetting, then load into a NucleoCassette.
7. Repeat step 6 for a second time, and calculate the average result.
8. Determine the total **Volume of Cell Suspension**.

**Note:** We use a 1,000  $\mu\text{L}$  pipette to establish this volume.

9. Calculate the number of viable cells in the cell suspension:

- ❖ Viable cells/mL =  $10 \times [(3 \times \text{Total Number of Cells/mL}) - (\text{Non-Viable Cells/mL})]$
- ❖ Total number of viable cells = Viable cells/mL  $\times$  Volume of Cell Suspension (mL)

**Worked example for a typical  $3.5 \times 10^6$  vial of Cytiva Plus Cardiomyocytes:**

Total Number of Cells/mL determined in step 6 =  $1.83 \times 10^5$

Non-Viable Cells/mL determined in step 5 =  $1.64 \times 10^5$

Viable cells/mL =  $10 \times [(3 \times 1.83 \times 10^5) - (1.64 \times 10^5)] = 3.85 \times 10^6$  cells/mL

Total number of viable cells =  $3.85 \times 10^6$  cells/mL  $\times$  1 mL =  $3.85 \times 10^6$  cells

**Note:** There are a total of  $5.5 \times 10^6$  cells in a  $3.5 \times 10^6$  vial of Cytiva Plus Cardiomyocytes. We routinely record a value of  $>70\%$  post thaw cell viability (i.e.  $3.85 \times 10^6$  viable cells).

## 5.5. Seeding Cytiva Plus Cardiomyocytes onto MEA plates (Time 3 hours)

Cells should be seeded at a density of  $6 \times 10^4$  viable cells in 4  $\mu\text{L}$  (i.e.  $1.5 \times 10^7$  viable cells/mL) over the recording electrode area.

**Note:** Assuming the volume of the cell suspension is  $\sim 1$  mL, and the number of viable cells is  $3.85 \times 10^6$ , the current concentration is  $\sim 3.85 \times 10^6$  viable cells/mL. Consequently, the next step involves centrifuging and re-suspending the cells from section 5.4 to achieve the appropriate cell concentration.

### Consumables

- RPMI 1640/B27 medium (filter sterilized)
- 50 mL tube
- 60 mL bottle

### Equipment

- Fibronectin coated well MEA plate
- 20  $\mu\text{L}$  pipette & sterile 1-20  $\mu\text{L}$  tips
- 200  $\mu\text{L}$  pipette & sterile 1-200  $\mu\text{L}$  tips
- 1,000  $\mu\text{L}$  pipette & sterile 1-1,000  $\mu\text{L}$  tips

### Protocol

Perform the following steps aseptically inside the biosafety cabinet (BSC).

1. Warm 25 mL RPMI 1640/B27 medium in a 60 mL bottle (sealed) in a  $37^\circ\text{C}$  water bath.
2. Centrifuge the cell suspension in a 50 mL tube at 300 g for 5 minutes at  $20^\circ\text{C}$ .
3. Carefully, using a 1,000  $\mu\text{L}$  pipette, remove the supernatant.

4. Re-suspend the cells in the residual liquid with gentle agitation.
5. Determine the total volume of the cell suspension.


**Note:** We use a 100  $\mu\text{L}$  pipette to establish this volume.

6. Dilute the cell suspension to  $1.5 \times 10^7$  viable cells/mL using warm RPMI 1640/B27 medium.

$$\text{Required volume (mL)} = \frac{\text{Viable cells}}{1.5 \times 10^7 \text{ viable cells/mL}}$$

**Note:** Assuming number of viable cells is  $3.85 \times 10^6$ , the final volume should be 260  $\mu\text{L}$ .

7. Remove fibronectin-coated well MEA plate from incubator.
8. Aspirate the fibronectin bead from each well of the MEA plate (use a pipette set to dispense 8  $\mu\text{L}$ ) and **immediately** replace with a 4  $\mu\text{L}$  bead of Cytiva Plus Cardiomyocyte suspension ( $1.5 \times 10^7$  viable cells/mL) over the recording electrode area (see Fig 2).




**Figure 2.** Fibronectin bead replaced with Cytiva Plus cardiomyocyte bead.

9. After seeding all wells, put the lid on the tissue culture dish and incubate the MEA plate in a standard cell culture incubator at  $37^\circ\text{C}$  for 2–3 hours. The sterile distilled water in the gaps between the wells will prevent the evaporation of the bead of cell suspension.
10. After 2–3 hours, gently add 500  $\mu\text{L}$  of warm RPMI 1640/B27 medium to the ‘corner’ of each well of the 12-well MEA plate or 150  $\mu\text{L}$  to the ‘corner’ of each well of the 48-well MEA plate. Pipette against the side of the well so as not to disturb the plated cells.
11. Gently add another 500  $\mu\text{L}$  of warm RPMI 1640/B27 medium to the ‘corner’ of each well of the 12-well MEA plate or 150  $\mu\text{L}$  to the ‘corner’ of each well of the 48-well MEA plate.

12. Gently add another 1,000  $\mu\text{L}$  of warm RPMI 1640/B27 medium to the 'corner' of each well of the 12-well MEA plate (i.e. total well volume is now 2 mL) or another 150  $\mu\text{L}$  to the 'corner' of each well of the 48-well MEA plate (i.e. total well volume is now 450  $\mu\text{L}$ ).
13. Aspirate the 500  $\mu\text{L}$  of sterile distilled  $\text{H}_2\text{O}$  from the gaps between the wells in the 12-well MEA plate or 200  $\mu\text{L}$  of sterile distilled  $\text{H}_2\text{O}$  from the gaps between the wells in the 48-well MEA plate.
14. Confirm cell attachment at this point by observing the MEA wells under the microscope if using a transparent plate (see Fig 3A and 3B).
15. Incubate in a standard cell culture incubator at  $37^\circ\text{C}$ , 5%  $\text{CO}_2$ .



**Figure 3A** – Wells containing Cytiva Plus Cardiomyocytes and media.



**Figure 3B** –Cytiva Plus Cardiomyocytes as viewed under a microscope, 5 days after seeding.

## 5.6. Media change on day 4 post-thaw (Time 20 minutes)

1. Warm 30 mL RPMI 1640/B27 medium in a 60 mL bottle in a  $37^\circ\text{C}$  water bath.
2. Perform the following steps aseptically inside the biosafety cabinet (BSC).
3. For the 12-well MEA plate, using a 1,000  $\mu\text{L}$  pipette, carefully aspirate 1,000  $\mu\text{L}$  of media from each well of the 12-well MEA plate, then gently add 1,000  $\mu\text{L}$  of warm RPMI 1640/B27 medium to the 'corner' of each well of the 12-well MEA plate, to give a total volume of 2 mL.
4. For the 48-well plate, using a 1,000  $\mu\text{L}$  pipette, carefully aspirate 250  $\mu\text{L}$  of media from each well of the 48-well MEA plate then using a 200  $\mu\text{L}$  pipette, gently add 200  $\mu\text{L}$  of warm RPMI 1640/B27 medium to the 'corner' of each well of the 48-well MEA plate, to give a total well volume of 400  $\mu\text{L}$ .
5. Perform MEA recordings 5–7 days after plating onto the MEA plate. User must determine the optimal culture time for this application.

## 6. Protocols for High Content Analysis (HCA)

### 6.1. Preparation of cell culture plates for HCA (Time 2½ hours)

#### Consumables

- Sterile distilled H<sub>2</sub>O
- Fibronectin (1 mg, human, BD Biosciences 354008) (at 4°C)
- 50 mL sterile tube
- D-PBS (sigma D8537)

#### Equipment

- 384-well cell culture plate (e.g. Greiner µClear 781091)\*
- 96-well cell culture plate (e.g. Greiner µClear 655090)\*
- 20 µL pipette & sterile 1–40 µL tips
- 1,000 µL pipette & sterile 1–1,000 µL tips

**\*Note:** Cytiva Plus Cardiomyocyte - ES-7002 (3.5x10<sup>6</sup> viable cells/vial) provides sufficient cardiomyocytes for **1 96-well plate (3.0x10<sup>4</sup> cells/well) or 1 384-well plate (9.0x10<sup>3</sup> cells/well).**

#### Protocol

Perform the following steps aseptically inside the biosafety cabinet (BSC).

1. Prepare a 1 mg/mL solution of fibronectin by adding 1 mL sterile distilled water to the 1 mg fibronectin.
2. Take 125 µL of this 1 mg/mL fibronectin solution and add it to 9,875 µL D-PBS in a 50 mL tube for a final concentration of 12.5 µg/mL fibronectin. Add 100 µL of 12.5 µg/mL fibronectin solution to each well of the 96-well cell culture plate, or 30 µL to each well of a 384-well cell culture plate.
3. Put the lid on the cell culture plate and incubate for 2 hours in a standard cell culture incubator at 37°C, 5% CO<sub>2</sub>.

### 6.2. Preparation of RPMI 1640/B27 medium for HCA (Time 20 minutes)

#### Consumables

- RPMI 1640 + Glutamine medium, 500 mL (at 4°C)
- B27 supplement, 10 mL (at -20°C)
- 10 mL sterile serological pipette

#### Equipment

- 500 mL 0.22 µm cellulose acetate filter unit (Corning 430769)
- Vacuum line
- Pipette gun

## Protocol

During the 2-hour fibronectin incubation step, prepare medium.

**Note:** Medium may be prepared immediately before thawing Cytiva Plus Cardiomyocytes or prepared and stored at 2–8°C and used within one week.

1. Thaw frozen 10 mL B27 supplement vial(s) for 10 minutes in a 37°C water bath. Do not incubate at 37°C or expose to light for extended periods of time. Perform the following steps aseptically inside the biosafety cabinet (BSC).
2. Wipe the required number of RPMI 1640 + Glutamine medium bottle(s) and B27 Supplement vial(s) with 70% isopropanol and transfer to a BSC.
3. Place the filtration unit into the BSC.
4. Using a 10 mL sterile serological pipette, add 10 mL of B27 supplement to 500 mL bottle of RPMI 1640 + Glutamine medium. Swirl bottle several times to mix.
5. Carefully pour 510 mL of the RPMI 1640/B27 medium into the reservoir of the filtration unit.
6. Place the lid on the filter unit.
7. Connect the filter unit to a vacuum source.
8. When filtration is complete, disconnect the vacuum source.
9. Detach the upper reservoir of the filtration unit.
10. Place sterile cap on the bottle portion of the filter unit.

**Note:** Store the medium at 2–8°C. Use within one week of preparation. Avoid repeated warming of the RPMI 1640/ B27 medium. Warm only the required volume of medium to complete the task.

## 6.3. Thawing Cytiva Plus Cardiomyocytes for HCA (Time 30 minutes)

### Consumables

- Cytiva Plus Cardiomyocytes- ES-7002
- RPMI 1640/B27 medium (filter sterilized)
- 50 mL sterile tube
- 10 mL sterile serological pipette

### Equipment

- Ice bucket with dry-ice
- Cryovial rack
- 37°C water bath
- Centrifuge
- 1000 µL pipette & sterile 1-1000 µL tips
- Pipette gun

## Protocol

During the 2-hour fibronectin incubation step, perform the following steps aseptically inside the biosafety cabinet (BSC).

1. Remove the cryovial from the cryostore and place onto dry-ice until ready to thaw.
2. Thaw cell suspension in a 37°C water bath with gentle agitation until ice crystals just disappear.

**Note:** Take care not to immerse the whole cryovial into the water bath. Avoid extended incubation at 37°C.

3. Wipe the outside of the cryovial with 70% isopropanol and transfer to BSC.
4. Carefully transfer the cell suspension into a sterile 50 mL centrifuge tube using a 1,000 µL pipette.
5. Rinse the inside of the cryovial with 1 mL of room temperature RPMI 1640/B27 and combine with the cell suspension drop-wise with gentle mixing.
6. Slowly (over the course of 2 minutes) add 8 mL of RPMI 1640/B27 to the 50 mL centrifuge tube.
7. Centrifuge at 300 g for 5 minutes at 20°C.
8. Carefully, using a 10 mL stripette, remove 8 mLs of the supernatant. Remove a further 1 mL with a 1,000 µL pipette taking care not to disturb the cell pellet. Re-suspend the cells in the residual liquid (approx. 1 mL) with gentle agitation.

#### 6.4. Determining Post-thaw Cell Viability of Cytiva Plus Cardiomyocytes for HCA (Time 30 minutes)

During the 2-hour fibronectin incubation step, determine the viable cell number and viable cell density using preferred method of choice. We use a NucleoCounter NC-100 (ChemoMetec) cell counter, the method for which is outlined below. However, this protocol could be adapted for other commercial cell counters.

##### Consumables

- RPMI 1640/B27 medium (filter sterilized)
- 4x1.5 mL tubes
- 4xNucleoCassettes™ (ChemoMetec)
- Reagent A100 (ChemoMetec)
- Reagent B (ChemoMetec)

##### Equipment

- 100 µL pipette & sterile 1-200 µL tips
- 1,000 µL pipette & sterile 1-1,000 µL tips
- NucleoCounter NC-100 (ChemoMetec)



## Protocol

Perform the following steps aseptically inside the biosafety cabinet (BSC).

1. Determine viable cell number using the NucleoCounter NC-100 (ChemoMetec). To do this, transfer 40  $\mu\text{L}$  of the cell suspension into a 1.5 mL tube. Add to the tube 360  $\mu\text{L}$  of RPMI 1640/B27 medium (i.e. now a 1:10 dilution). Gently agitate to achieve an even cell suspension.
2. The remaining steps can be performed outside of the BSC.
3. Transfer 100  $\mu\text{L}$  of 1:10 diluted cell suspension into each of three more 1.5 mL tubes.
4. Use a 100  $\mu\text{L}$  of 1:10 diluted cell suspension sample to calculate the number of **Non-Viable Cells/mL** (i.e. directly load the sample into a NucleoCassette).
5. Repeat step 4 for a second time, and calculate the average result.
6. Use another 100  $\mu\text{L}$  of 1:10 diluted cell suspension sample to calculate the **Total Number of Cells/mL**. To do this, add 100  $\mu\text{L}$  of Reagent A100, and then 100  $\mu\text{L}$  of Reagent B, to the 100  $\mu\text{L}$  of 1:10 diluted cell suspension sample. Mix by pipetting, then load into a NucleoCassette.
7. Repeat step 6 for a second time, and calculate the average result.
8. Determine the total **Volume of Cell Suspension**.

**Note:** We use a 1,000  $\mu\text{L}$  pipette to establish this volume.

9. Calculate the number of viable cells in the cell suspension:
  - ❖  $\text{Viable cells/mL} = 10 \times [(3 \times \text{Total Number of Cells/mL}) - (\text{Non-Viable Cells/mL})]$
  - ❖  $\text{Total number of viable cells} = \text{Viable cells/mL} \times \text{Volume of Cell Suspension (mL)}$

### Worked example for a typical $3.5 \times 10^6$ vial of Cytiva Plus Cardiomyocytes:

**Total Number of Cells/mL** determined in step 6 =  $1.83 \times 10^5$

**Non-Viable Cells/mL** determined in step 5 =  $1.64 \times 10^5$

$\text{Viable cells/mL} = 10 \times [(3 \times 1.83 \times 10^5) - (1.64 \times 10^5)] = 3.85 \times 10^6 \text{ cells/mL}$

$\text{Total number of viable cells} = 3.85 \times 10^6 \text{ cells/mL} \times 1 \text{ mL} = 3.85 \times 10^6 \text{ cells}$

**Note:** There are a total of  $5.5 \times 10^6$  cells in a  $3.5 \times 10^6$  vial of Cytiva Plus Cardiomyocytes and we routinely record a value of >70% post thaw cell viability (i.e.  $3.85 \times 10^6$  viable cells).

## 6.5. Seeding Cytiva Plus Cardiomyocytes into cell culture plates for HCA (Time 30 minutes)

Cells should be seeded at a density of  $3.6 \times 10^4$  viable cells in 200  $\mu$ L RPMI 1640/B27 per well of a 96-well cell culture plate, or  $9 \times 10^3$  viable cells in 50  $\mu$ L RPMI 1640/B27 per well of a 384-well cell culture plate (i.e.  $1.8 \times 10^5$  viable cells/mL).

**Note:** Assuming the volume of the cell suspension is  $\sim 1$  mL, and the number of viable cells is  $3.85 \times 10^6$ , the current concentration is  $\sim 3.85 \times 10^6$  viable cells/mL. Consequently, the next step involves diluting the cell suspension to achieve the appropriate cell concentration.

### Consumables

- RPMI 1640/B27 medium (filter sterilized)
- 60 mL bottle

### Equipment

- Fibronectin coated cell culture plate
- 200  $\mu$ L pipette & sterile 1-200  $\mu$ L tip

### Protocol

Perform the following steps aseptically inside the biosafety cabinet (BSC).

1. Warm 25 mL RPMI 1640/B27 medium in a sterile 60 mL bottle (sealed) in a 37°C water bath.
2. Dilute the cell suspension to  $1.8 \times 10^5$  viable cells/mL using warm RPMI 1640/B27 medium.


$$\text{Required volume (mL)} = \frac{\text{Viable cells}}{1.8 \times 10^5 \text{ viable cells/mL}}$$

**Example:** A total viable count of  $3.85 \times 10^6$  cells would be divided by  $1.8 \times 10^5$  viable cells/mL to obtain a final required volume of 21.4 mL.

3. Remove fibronectin-coated cell culture plate from the incubator.
4. Aspirate the fibronectin from each well of the cell culture and **immediately** replace with 200  $\mu$ L of Cytiva Plus Cardiomyocyte suspension (i.e.  $1.8 \times 10^5$  viable cells/mL) per well of a 96-well culture plate, or 50  $\mu$ L per well of a 384-well culture plate.
5. After seeding all the required wells, put the lid on the cell culture plate and incubate in a standard cell culture incubator at 37°C.

## 6.6. Media change on day 4 post-thaw (Time 20 minutes)

1. Warm 30 ml RPMI 1640/B27 medium in a 60 ml bottle (sealed) in a 37°C water bath.
2. Perform the following steps aseptically inside the biosafety cabinet (BSC).
3. After 96 hours, flick the media out of the plate onto clean tissues and replace with 40 µL of warmed RPMI 1640/B27 medium into each well so as not to disturb the seeded cells. Incubate plate at 37°C, 5% CO<sub>2</sub>.
4. Repeat a media change at day 7 replacing the seeded medium with fresh pre-warmed RPMI 1640/B27 medium.
5. Perform HCA 7–8 days after seeding in cell culture plates. It is the responsibility of the user to determine the optimal culture time for this application.



**Figure 4** – Bright field image of Cytiva Plus Cardiomyocytes seeded in microplates for HCA at 8 days post thaw.

## 7. Protocols for Manual Patch Clamp (MPC)

### 7.1. Sterilization of glass coverslips (Time 4 hours)

#### Consumables

- Borosilicate glass coverslips (13 mm)
- Isopropanol
- 24-well cell culture plate
- 50 mL sterile tube

#### Equipment

- Sterile forceps

#### Protocol

Perform the following steps aseptically inside the biosafety cabinet (BSC).

1. Sterilize the coverslips by soaking in isopropanol for a minimum of 2 hours in a sealed sterile 50 mL tube. Ensure complete wetting of both sides of each coverslip.
2. Using sterile forceps, remove a coverslip from the isopropanol and shake the excess isopropanol from the coverslip. Place the coverslip into the well of a sterile 24-well plastic cell culture plate.
3. Repeat steps 1–3 for the number of wells required.
4. Let the coverslips air-dry in a BSC for a minimum of 2 hours.

### 7.2. Preparation of 1:2 diluted Matrigel aliquots (10 minutes with O/N thaw)

#### Consumables

- Matrigel (BD Biosciences 356231, -20°C)
- KO-DMEM (2–8°C)
- 10 mL sterile serological pipette (2–8°C)
- 10x50 mL sterile tubes (2–8°C)

#### Equipment

- Pipette gun

#### Protocol

Do not allow the Matrigel solution to reach room temperature. Keep the solution, pipettes, and KO-DMEM cold at all stages of handling. Avoid repeated freeze-thawing of diluted Matrigel aliquots. Perform the following steps aseptically inside the BSC.

1. Slowly thaw Matrigel at 2–8°C **overnight** to avoid the formation of a gel.
2. Cool a sterile 10 mL serological pipette by drawing and releasing 10 mL of cold KO-DMEM into the pipette repeatedly without removing the pipette from the bottle of KO-DMEM.
3. Add 10 mL of cold KO-DMEM to the vial containing 10 mL Matrigel.
4. Working quickly, mix the Matrigel and KO-DMEM with a 10 mL pipette, avoiding the formation of bubbles.
5. Aliquot 2 mL of diluted Matrigel into pre-chilled sterile 50 mL tube; store at -20°C until required. Diluted Matrigel solution is stable for 3 months when stored at -20°C.

### 7.3. Preparation of 1:30 diluted Matrigel (Time 10 minutes)

#### Consumables

- 2 mL 1:2 diluted Matrigel (-20°C)
- KO-DMEM (2–8°C)
- 5 mL sterile serological pipette (2–8°C)
- 25 mL sterile serological pipette (2–8°C)

#### Equipment

- Pipette gun

#### Protocol

Do not allow the Matrigel solution to reach room temperature. Keep the solution and pipettes cold at all stages of handling. Perform the following steps aseptically inside the biosafety cabinet (BSC).

1. Slowly thaw a 2 mL 1:2 diluted Matrigel aliquot prepared in step 7.2 at 4°C for at least 2 hours to avoid the formation of a gel. Once thawed, transfer to the BSC.
2. Cool a sterile 5 mL pipette by drawing 5 mL of cold KO-DMEM into the pipette.
3. Dilute the 2 mL 1:2 diluted Matrigel aliquot with 5 mL cold KO-DMEM. Carefully mix Matrigel solution, avoiding the formation of bubbles.
4. Cool a sterile 25 mL pipette by drawing 25 mL of cold KO-DMEM into the pipette.
5. Add a further 23 mL cold KO-DMEM (for a final dilution of 1:30). Carefully mix Matrigel solution, avoiding the formation of bubbles.

### 7.4. Preparation of Matrigel-coated coverslips (Time 10 minutes)

#### Consumables

- 30 mL 1:30 diluted Matrigel (4°C)
- 24-well cell culture plate containing glass coverslips
- KO-DMEM (2–8°C)
- 10 mL sterile serological pipette (2–8°C)
- 10 × 50 mL tubes

#### Equipment

- Pipette gun

#### Protocol

Do not allow the Matrigel solution to reach room temperature. Keep the solution, pipettes and KO-DMEM cold at all stages. Avoid repeated freeze-thawing of diluted Matrigel aliquots. Perform the following steps aseptically inside the biosafety cabinet (BSC).

1. Add 400 µL of 1:30 diluted Matrigel to each well of the 24-well cell culture plate containing a glass coverslip.
2. Incubate the cell culture plate overnight at 2–8°C before use.
3. Matrigel coated vessels are stable for 10 days when stored at 2–8°C.

## 7.5. Preparation of RPMI 1640/B27 medium for MPC (Time 20 minutes)

### Consumables

- RPMI 1640 + Glutamine medium, 500 mL (at 4°C)
- B27 supplement, 10 mL (at 4°C)
- 10 mL sterile serological pipette

### Equipment

- 500 mL 0.22 µm cellulose acetate filter unit (Corning 430769)
- Vacuum line
- Pipette gun

### Protocol

Perform the following steps aseptically inside the biosafety cabinet (BSC).

1. Using a 10 mL sterile serological pipette, add 10 mL of B27 supplement to 500 mL bottle of RPMI 1640 + Glutamine medium. Swirl bottle several times to mix.
2. Place the filtration unit into the BSC.
3. Carefully pour 510 mL of the RPMI/B27 medium into the reservoir of the filtration unit.
4. Place the lid on the filter unit.
5. Connect the filter unit to a vacuum source.
6. When filtration is complete, disconnect the vacuum source.
7. Detach the upper reservoir of the filtration unit.
8. Place sterile cap on the bottle portion of the filter unit.

**Note:** Store the medium at 2–8°C. Use within one week of preparation. Avoid repeated warming of the RPMI 1640/B27 medium. Warm only the required volume of medium to complete the task.

## 7.6. Thawing Cytiva Plus Cardiomyocytes for MPC (Time 30 minutes)

### Consumables

- 24-well cell culture plate containing Matrigel-coated coverslips (2–8°C)
- Cytiva Plus Cardiomyocytes (ES-7001 -  $1.0 \times 10^6$  cells/vial)
- RPMI 1640/B27 medium (filter sterilized)
- 50 mL sterile tube
- 10 mL sterile serological pipette

### Equipment

- Ice bucket with dry-ice
- Cryovial rack
- 37°C water bath
- Centrifuge
- 1000  $\mu$ L pipette & sterile 1–1000  $\mu$ L tips
- Pipette gun

### Protocol

Perform the following steps aseptically inside the biosafety cabinet (BSC).

1. Before thawing the cryovial of cells, take the 24-well cell culture plate containing the Matrigel-coated coverslips out of 2–8°C storage. Leave the 24-well cell culture plate at room temperature for 1 hour before seeding the cells.
2. Remove the cryovial from the cryostore and place onto dry-ice until ready to thaw.
3. Thaw cell suspension in a 37°C water bath with gentle agitation until ice crystals just disappear.

**Note:** Take care not to immerse the whole cryovial into the water bath. Avoid extended incubation at 37°C.

4. Wipe the outside of the cryovial with 70% isopropanol and transfer to BSC.
5. Carefully transfer the cell suspension into a sterile 50 mL centrifuge tube using a 1,000  $\mu$ L pipette.
6. Rinse the inside of the cryovial with 1 mL of room temperature RPMI 1640/B27 and combine with the cell suspension drop-wise with gentle mixing.
7. Slowly (over 2 minutes) add 8 mL of RPMI 1640/B27 to a 50 mL centrifuge tube.
8. Centrifuge at 300 g for 5 minutes at 20°C.
9. Carefully, using a 10 mL stripette remove 8 mLs of the supernatant. Remove a further 1 mL with a 1,000  $\mu$ L pipette taking care not to disturb the cell pellet. Re-suspend the cells in the residual liquid (approx. 1 mL) with gentle agitation.

## 7.7. Determining Post-thaw Cell Viability of Cytiva Plus Cardiomyocytes for MPC (Time 30 minutes)

Determine the viable cell number and viable cell density using preferred method of choice. We use a NucleoCounter NC-100 (ChemoMetec) cell counter, the method for which is outlined below. However, this protocol could be adapted for other commercial cell counters.

### Consumables

- RPMI 1640/B27 medium (filter sterilized)
- 4×1.5 mL tubes
- 4×NucleoCassettes™ (ChemoMetec)
- Reagent A100 (ChemoMetec)
- Reagent B (ChemoMetec)

### Equipment

- 100 µL pipette & sterile 1-200 µL tips
- 1,000 µL pipette & sterile 1-1,000 µL tips
- NucleoCounter NC-100 (ChemoMetec)

### Protocol

Perform the following step aseptically inside the biosafety cabinet (BSC).

1. Determine viable cell number using the NucleoCounter NC-100 (ChemoMetec). To do this, transfer 40 µL of the cell suspension into a 1.5 mL tube. Add to the tube 360 µL of RPMI 1640/B27 medium (i.e. now a 1:10 dilution). Gently agitate to achieve an even cell suspension.
2. The remaining steps can be performed outside of the BSC.
3. Transfer 100 µL of 1:10 diluted cell suspension into each of three more 1.5 mL tubes.
4. Use a 100 µL of 1:10 diluted cell suspension sample to calculate the number of **Non-Viable Cells/mL** (i.e. directly load the sample into a NucleoCassette).
5. Repeat step 4 for a second time, and calculate the average result.
6. Use another 1×100 µL of 1:10 diluted cell suspension sample to calculate the **Total Number of Cells/mL**. To do this, add 100 µL of Reagent A100, and then 100 µL of Reagent B, to the 100 µL of 1:10 diluted cell suspension sample. Mix by pipetting, then load into a NucleoCassette.
7. Repeat step 6 for a second time, and calculate the average result.
8. Determine the total **Volume of Cell Suspension**.



**Note:** We use a 1,000  $\mu\text{L}$  pipette to establish this volume.

9. Calculate the number of viable cells in the cell suspension:

- ❖ Viable cells/mL =  $10 \times [(3 \times \text{Total Number of Cells/mL}) - (\text{Non-Viable Cells/mL})]$
- ❖ Total number of viable cells = Viable cells/mL  $\times$  Volume of Cell Suspension (mL)

**Worked example for a typical  $1.0 \times 10^6$  vial of Cytiva Plus Cardiomyocytes:**

Total Number of Cells/mL determined in step 6 =  $5.33 \times 10^4$

Non-Viable Cells/mL determined in step 5 =  $4.80 \times 10^4$

Viable cells/mL =  $10 \times [(3 \times 5.33 \times 10^4) - (4.80 \times 10^4)] = 1.12 \times 10^6$  cells/mL

Total number of viable cells =  $1.12 \times 10^6$  cells/mL  $\times$  1 mL =  $1.12 \times 10^6$  cells

**Note:** There are a total of  $1.6 \times 10^6$  cells in a  $1.0 \times 10^6$  vial of Cytiva Plus Cardiomyocytes. We routinely record a value of  $>70\%$  post thaw cell viability (i.e.  $1.12 \times 10^6$  viable cells).

## 7.8. Seeding Cytiva Plus Cardiomyocytes onto coverslips (Time 30 minutes)

Cells should be seeded at a density of  $4.6 \times 10^4$  viable cells in 650  $\mu\text{L}$  RPMI 1640/B27 per well of a 24-well cell culture plate (i.e.  $7.1 \times 10^4$  viable cells/mL).

**Note:** Assuming the volume of the cell suspension is  $\sim 1$  mL, and the number of viable cells is  $1.12 \times 10^6$ , the concentration is  $\sim 1.12 \times 10^6$  viable cells/mL. Consequently, the next step involves diluting the cell suspension to achieve the appropriate cell concentration.

### Consumables

- RPMI 1640/ medium (filter sterilized)
- 60 mL bottle

### Equipment

- 24-well culture plate containing Matrigel-coated coverslips
- 200  $\mu\text{L}$  pipette & sterile 1-200  $\mu\text{L}$  tips
- $37^\circ$  water bath

### Protocol

Perform the following steps aseptically inside the biosafety cabinet (BSC).

1. Warm 25 mL RPMI 1640/B27 medium in a sterile 60 mL bottle (sealed) in a  $37^\circ\text{C}$  water bath.
2. Dilute the cell suspension to  $7.1 \times 10^4$  viable cells/mL using warm RPMI 1640/B27 medium.


$$\text{Required volume (mL)} = \frac{\text{Viable cells}}{7.1 \times 10^4 \text{ viable cells/mL}}$$

**Example:** A total viable count of  $1.12 \times 10^6$  cells would be divided by  $7.1 \times 10^4$  viable cells/mL to obtain a final required volume of 15.8 mL.

3. Aspirate the Matrigel from each well of the 24-well cell culture plate and **immediately** replace with 650  $\mu\text{L}$  of Cytiva Plus Cardiomyocyte suspension (i.e.  $7.1 \times 10^4$  viable cells/mL).
4. After seeding all the required wells, put the lid on the cell culture plate and incubate in a standard cell culture incubator at  $37^\circ\text{C}$ , 5%  $\text{CO}_2$ .

### 7.9. Media change on day 4 post-thaw (Time 20 minutes)

1. Warm 30 mL RPMI 1640/B27 medium in a 60 mL bottle (sealed) in a  $37^\circ\text{C}$  water bath.
2. Perform the following steps aseptically inside the biosafety cabinet (BSC).
3. After 96 hours, carefully aspirate off 300  $\mu\text{L}$  of the medium from each well, leaving 350  $\mu\text{L}$  residual medium and replace with 300  $\mu\text{L}$  of warmed RPMI 1640/B27 medium so as not to disturb the seeded cells. Incubate plate at  $37^\circ\text{C}$ , 5%  $\text{CO}_2$ .
4. Repeat a media change at day 7 replacing half the seeded medium with fresh pre-warmed RPMI 1640/B27 medium.
5. Perform manual patch clamp 4–7 days after seeding in cell culture plates. It is the responsibility of the user to determine the optimal culture time for this application.



**Figure 5** – Examples of cell morphologies patched manually.

## 8. Troubleshooting

### Problems

### Solutions

---

Poor cardiomyocyte viability

Check storage/shipping conditions. Follow pack leaflet to revive cells from cryostore. Follow recommended procedure for thawing and dilution of cardiomyocytes as improper handling of Cytiva Plus Cardiomyocytes could cause low viability. Ensure pre-warmed medium is added drop-wise to the cardiomyocytes.

Poor cell attachment

Make sure plates are tissue culture treated and sterile. Follow recommended procedure for preparation and storage of fibronectin plates. Aspirate fibronectin just before seeding. Be careful not to touch the bottom of plate. Make sure the surface of the well has not been marked during the removal of fibronectin. Do not let the plate dry out.

Place plate into incubator on a flat surface for 2 hours.

Make sure viable count number is used to seed.

Follow feeding steps as outlined.

Ensure RPMI 1640/B27 is: a) warmed before use; b) used within 7 days of preparation; c) not exposed to light for extended periods; d) used within 7 days of preparation, and; e) not warmed repeatedly.

Ensure that Matrigel coverage is uniform by microscopic inspection. Matrigel coated plates should be warmed to room temperature for at least 30 minutes before use.

## Problems

## Solutions

Uneven coating of Matrigel

Ensure culture vessels are stored on a level surface. Use only cell culture treated plastic ware.

Aggregation of Matrigel

Follow the recommended procedure for the preparation of Matrigel coating of plates.

Avoid warming of Matrigel. Keep Matrigel on ice during processing if necessary.

Avoid room temperature plastic coming in contact with Matrigel.

Poor ion channel activity

Allow at least 72-hour post thaw before use depending on the application. Ion channel activity will increase as cells recover from thaw.

Ensure correct seeding density of Cytiva Plus for the application.

Poor performance on MEA

Do not allow the droplets of cell suspension (5.5.9) to dry out. The cell droplet should still be visible post 2–3 hour incubation prior to adding medium to the well (5.5.10).


Take care when adding the medium to the well so as not to disturb the seeded cells (5.5.10).

If using re-usable MEA plates, ensure that the MEA plates are:

(A) Sterile and the plating surface is dry and free of cellular debris.

(B) Maintained and stored as directed by the manufacturer.

## 9. Product Protocol Flowchart



**ASCENDANCE™**